



POWDERED

POWDER-FREE  
EXAM GRADE

TEXTURED  
GRIP

SINGLE  
USE ONLY

EXTENDED  
LENGTH

## Ultra Thick Disposable Gloves

- Powdered latex or powder-free exam grade latex
- Thickness: 14 mil
- 12" length for added protection
- Beaded cuff for easy donning
- Fully textured for superior grip
- Strong yet soft and flexible
- Single use only

### Powder-Free

Item No.	Size	Retail Qty	Case Qty	UPC	ITF10
6602-20	Medium	50 Gloves	10 Boxes	781311660229	10781311660226
6603-20	Large	50 Gloves	10 Boxes	781311660328	10781311660325
6604-20	X-Large	50 Gloves	10 Boxes	781311660427	10781311660424
6605-20	2X-Large	50 Gloves	10 Boxes	781311660526	10781311660523

### Powdered

Item No.	Size	Retail Qty	Case Qty	UPC	ITF10
6602	Medium	50 Gloves	10 Boxes	781311066021	10781311066028
6603	Large	50 Gloves	10 Boxes	781311066038	10781311066035
6604	X-Large	50 Gloves	10 Boxes	781311066045	10781311066042
6605	2X-Large	50 Gloves	10 Boxes	781311066052	10781311066059

Item No.	Size	Retail Qty	Case Qty	UPC	ITF10
6603-22	Large	2 Pairs	25 Packs	781311663220	10781311663227
6604-22	X-Large	2 Pairs	25 Packs	781311664227	10781311664224



## Palm Width (mm)

Size	Palm Width (mm)
Medium	97 ± 5
Large	107 ± 5
X-Large	117 ± 5
2X-Large	125 ± 5

## Length / Thickness (single wall)

Dimension	Value
Length (mm)	300 ± 5
Cuff Thickness (mm)	0.20 ± 0.02
Palm Thickness (mm)	0.32 ± 0.02
Finger Thickness (mm)	0.36 ± 0.02

## Tensile Strength

Before Aging (MPa)	min. 22
After Aging (MPa)	min. 18

## Elongation

Before Aging (%)	750
After Aging (%)	680